













EasyShock 2.0 Transmitter energy density list(mJ/mm2)

Energy (mJ)	8mm (1#)	15mmA (2#)	15mmB (3#)	20mm (4#)	36mm (5#)	ED1 (6#)	ED2 (7#)
Treatment head picture							
210	4.18	1.11	1.17	0.65	0.21	2.68	2.86
200	3.98	1.05	1.11	0.62	0.20	2.55	2.72
190	3.78	1.00	1.06	0.59	0.19	2.42	2.59
180	3.58	0.95	1.00	0.56	0.18	2.29	2.45
170	3.38	0.89	0.94	0.53	0.17	2.17	2.31
160	3.18	0.84	0.89	0.50	0.16	2.04	2.18
150	2.99	0.79	0.83	0.47	0.15	1.91	2.04
140	2.79	0.74	0.78	0.44	0.14	1.78	1.90
130	2.59	0.68	0.72	0.40	0.13	1.66	1.77
120	2.39	0.63	0.67	0.37	0.12	1.53	1.63
110	2.19	0.58	0.61	0.34	0.11	1.40	1.50
100	1.99	0.53	0.56	0.31	0.10	1.27	1.36
90	1.79	0.47	0.50	0.28	0.09	1.15	1.22
80	1.59	0.42	0.44	0.25	0.08	1.02	1.09
70	1.39	0.37	0.39	0.22	0.07	0.89	0.95
60	1.19	0.32	0.33	0.19	0.06	0.76	0.82
50	1.00	0.26	0.28	0.16	0.05	0.64	0.68
40	0.80	0.21	0.22	0.12	0.04	0.51	0.54
30	0.60	0.16	0.17	0.09	0.03	0.38	0.41

Shock wave therapy transmitter specifications

Transmitter No	Pictures	Size	Inddications	Depth
1#		8mm	Acupuncture therapy, thendon disorders, calcified shoulder,heel spur	0-60mm
2#		15mm(concave)	Soft tissue repair,Deep seated location , Calcifications	0-80mm
3#		15mm(Convex)	Pain therapy in tendon disorders,calcified shoulder, heel spur,and specific muscle points	0-50mm
4#		20mm	Myofascial trigger therapy, muscle and connective tissue, back muscles, dissolution of fibrosis, treatment of deep-seated muscles	0-80mm
5#		36mm	Myofascial trigger therapy, muscle and connective tissue, back muscles, skin elasticity, trigger-diagnose,dissolution of fibrosis, treatment of deep-seated muscles	0-50mm